


MS Spectrum



**Sample Information Item**  
 Date and Time : 2005-4-19 19:23:03  
 User : FJS  
 Sample : RP10452(GL-11)  
 Inj. Volume : 1.000  
 MW : 1053.17

**Probe:** ESI  
 Nebulizer Gas Flow: 1.5L/min  
 CDL: -20.0v  
 CDL Temp: 250°C  
 Block Temp: 200°C

**Probe bias:** +4.5kv  
 Detector: 1.5kv  
 T. Flow: 0.2ml/min  
 B. conc: 50%H2O/50%ACN