


MS Spectrum



Sample Information Item

Date and Time : 2005-1-17 18:32:13
User : FJS
Sample : RP10828(CY-10-1)
Inj. Volume : 1.000
MW : 1248.32

Cat. No.: RP10828

Lot No.: PJ00129504

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H₂O/50%ACN