


Positive

M/S Spectrum



Sample Information Item

Date and Time : 2005-3-4 14:29:57
 User : xie
 Sample : RP11291 SY-40
 Inj. Volume : 1.000
 MW : 4535.17

Cat. No.: RP11291

Lot No.: PJ00209510

Probe: ESI
 Nebulizer Gas Flow: 1.5l/min
 CDL: -20.0v
 CDL Temp: 250°C
 Block Temp: 200°C

Probe bias: +4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. conc: 50%H2O/50%ACN