


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-1-14 16:21:00  
User : xie  
Sample : RP11604(AA-11)  
Inj. Volume : 1.000  
MW : 1313.58

Cat. No.: RP11604  
Lot No.: PJ00521504

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp: 250°C  
Block Temp: 200°C

Probe bias: +4.5kv  
Detector: 1.5kv  
T. Flow: 0.2ml/min  
B. conc: 50%H2O/50%ACN