


MS Spectrum



Sample Information Item

Date and Time : 2005-1-26 12:21:40
User : FJS
Sample : RP13279(AE-6)
Inj. Volume : 1.000
MW : 673.79

Cat. No.: RP13279
Lot No.: PL00127501

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H₂O/50%ACN