


Sample Description:
Summary:

Product Name: RP10061 (CG-9-NH2)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
 A B
0.01min 7% 93%
25.0min 35% 65%
25.1min 100% 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10061
Lot No.: PJ00913505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.217	20091.588	171271.188	3.7143
2		16.522	371002.938	4439891.500	96.2857
Total			391094.525	4611162.688	100.0000