


Sample Description:
Summary:

Product Name : RP10094 (NV-41-NH2)
Column : 4.6*250mm, Hypersil-ODS2 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient

Cat. No.: RP10094
Lot No.: PJ00813505

	A	B
0.0min	38%	62%
25.0min	52%	48%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.722	2901.871	103701.313	0.8537
2		13.152	8484.676	251910.484	2.0738
3		14.030	301139.750	11585851.000	95.3781
4		14.030	13748.821	182073.781	1.4989
5		17.003	1197.427	23748.711	0.1955
Total			327472.545	12147285.289	100.0000