


Sample Description:  
Summary:

Product Name: RP10100 (FI-38-NH2)  
Column : 4.6×250mm, Sinochrom ODS-BP 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
0.01min      A      B  
                 35%      65%  
25.0min      40%      60%  
25.1min      100%      0%  
30.0min           STOP  
Flow rate : 1.0ml/min  
Wavelength : 220nm

Cat. No.: RP10100  
Lot No.: PJ00713505



**Results**

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.873	891.392	10050.854	0.2117
2		15.325	1289.079	42121.410	0.8873
3		16.985	1105.270	37915.367	0.7987
4		18.152	90149.305	4549613.500	95.8349
5		18.152	8.464	107644.781	2.2675
Total			93443.510	4747345.912	100.0000