


Sample Description:

Product Name : RP10101 (FI-38-NH2)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10101
 Lot No.: PJ00106506

	A	B
0.01min	35%	65%
25.0min	43%	57%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		8.730	343.265	10065.702	0.3230
2		9.807	675.373	18106.598	0.5811
3		12.095	811.933	34940.605	1.1213
4		13.985	87125.930	2985975.000	95.8221
5		13.985	3694.511	67075.875	2.1525
Total			92651.012	3116163.780	100.0000