


Sample Description:

Product Name : RP10104 (DV-40-NH2)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10104
 Lot No.: PJ00101506

	A	B
0.01min	40%	60%
25.0min	45%	55%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		9.400	318.178	2460.300	0.0417
2		11.150	2848.751	30527.496	0.5177
3		11.713	10145.607	180789.891	3.0658
4		12.713	781.228	9404.604	0.1595
5		15.748	109082.523	5606898.000	95.0821
6		15.748	2915.069	60199.090	1.0209
7		19.298	461.713	6624.806	0.1123
Total			126553.069	5896904.186	100.0000