


Sample Description:  
Summary:

Product Name : RP10178(RM-11-NH2-1)  
Column : 4.6\*250mm, Kromasil C18 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
0.01min A 23% B 77%  
25.0min A 40% B 60%  
25.1min A 100% B 0%  
30.0min STOP  
Flow rate : 1.0ml/min  
Wavelength : 220nm

Cat. No.: RP10178  
Lot No.: PJ00113505



**Results**

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		8.937	834.530	13639.600	0.3713
2		15.902	497.813	7116.500	0.1937
3		17.057	5607.325	87861.680	2.3917
4		17.530	120728.758	3565009.000	97.0433
Total			127668.426	3673626.780	100.0000