

Sample Description:


Product Name : RP10207 (RNL-11-NH2)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10207
 Lot No.: PJ01025504

	A	B
0.01min	27%	73%
25.0min	35%	65%
25.1min	100%	0%

STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		15.838	2214.913	32767.697	0.7937
2		16.347	207947.484	3958857.000	95.8943
3		17.063	491.536	5197.598	0.1259
4		17.918	6810.631	100889.898	2.4438
5		18.510	555.048	5260.202	0.1274
6		22.050	1578.052	25382.699	0.6148
Total			219597.665	4128355.095	100.0000