


Sample Description:

Product Name : RP10208(RM-11-NH2)  
 Column : 4.6\*250mm, Kromasil C18 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient  
     0.01min      A                      B  
     25.0min     20%                    80%  
     25.1min     45%                    55%  
     25.1min     100%                  0%  
     30.0min                            STOP  
 Flow rate : 1.0ml/min  
 Wavelength : 220nm

Cat. No.: RP10208  
 Lot No.: PJ00729504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		9.278	1380.397	11340.446	0.2021
2		17.197	2260.861	19498.852	0.3474
3		17.587	361206.781	5563307.500	99.1322
4		18.380	1557.178	17863.613	0.3183
Total			366405.218	5612010.411	100.0000