


Sample Description:
Summary:

Product Name : RP10214(TI-40)
Column : 4.6*250mm, Hypersil-ODS2 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient

Cat. No.: RP10214
Lot No.: PJ00415506

	A	B
0.01min	33%	67%
25.0min	60%	40%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		12.995	350.276	5730.012	0.2368
2		15.088	701.289	7313.160	0.3023
3		15.532	3199.552	61461.961	2.5403
4		15.978	126964.125	2344946.000	96.9206
Total			131215.242	2419451.133	100.0000