


Sample Description:

Product Name : RP10227[SC-25(0)]
 Column : 4.6*250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10227
 Lot No.: PJ00126507

	A	B
0.01min	25%	75%
25.0min	33%	67%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.567	7260.455	127711.680	0.3532
2		17.520	12356.036	353995.625	0.9790
3		17.950	44648.773	678081.500	1.8753
4		18.482	665736.063	34998580.000	96.7925
Total			730001.327	36158368.805	100.0000