


Sample Description:

Product Name : RP10229(SE-12)
 Column : 4.6*250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 8% 92%
 25.0min 30% 70%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP10229
 Lot No.: PJ00825504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.342	20459.139	202974.750	1.2207
2		13.047	4113.803	56213.816	0.3381
3		13.195	5676.472	58526.270	0.3520
4		13.907	3290.071	23731.621	0.1427
5		14.063	17954.934	185718.875	1.1169
6		14.552	871316.063	16100644.000	96.8296
Total			922810.480	16627809.332	100.0000