


Sample Description:
Summary:

Product Name : RP10315(RG-13)
Column : 4.6*250mm, Kromasil C18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A 10% B 90%
25.0min A 35% B 65%
25.1min A 100% B 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10315
Lot No.: PJ00513505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		18.418	683.365	8895.651	0.4399
2		18.837	3088.126	30875.705	1.5268
3		19.155	112232.898	1982478.000	98.0333
Total			116004.389	2022249.356	100.0000