


Sample Description:
Summary:

Product Name: RP10339(HK-38-NH2)
Column : 4.6×250mm, Kromasil c18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water

Cat. No.: RP10339
Lot No.: PJ00613505

Gradient
 A B
0.01min 20% 80%
25.0min 35% 65%
25.1min 100% 0%
30.0min STOP

Flow rate : 1.0ml/min
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		8.940	4927.901	46787.902	0.5013
2		17.697	13102.563	180916.469	1.9384
3		18.363	250807.875	9105591.000	97.5603
Total			268838.339	9333295.371	100.0000