

Sample Description:
Summary:


Product Name: RP10342(DA-29-2)
 Column : 4.6×250mm, Kromasil c18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water

Cat. No. : RP10342
 Lot No. : PX00125502

Gradient

	A	B
0.01min	24%	76%
25.0min	45%	55%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.:	Peak ID	Ret Time	Height	Area	Conc.
1		13.587	2307.869	18863.301	0.2707
2		13.912	1394.888	14401.555	0.2067
3		16.260	1085.182	8868.703	0.1273
4		17.237	1381.544	17246.297	0.2475
5		18.237	329032.531	6742150.000	96.7610
6		20.690	1575.565	30784.055	0.4418
7		22.412	7482.971	124056.992	1.7804
8		24.208	1087.711	11467.255	0.1646
Total			345348.260	6967838.157	100.0000