


Sample Description:
Summary:

Product Name: RP10399 (AY-36-2)
Column : 4.6×250mm, Kromasil C18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A B
 25% 75%
25.0min 45% 55%
25.1min 100% 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10399
Lot No.: PJ00213506



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		10.500	1972.992	22706.592	0.3258
2		13.658	3790.678	39235.551	0.5630
3		14.602	544162.875	6843255.500	98.2007
4		14.602	4685.169	33245.094	0.4771
5		15.420	2815.305	30202.496	0.4334
Total			557427.019	6968645.232	100.0000