


Sample Description:

Product Name : RP10415(RM-28-NH2)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10415
 Lot No.: PJ00108506

	A	B
0.01min	25%	75%
25.0min	50%	50%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.287	1651.891	14119.949	0.3994
2		12.438	1019.369	8090.290	0.2288
3		12.798	261777.297	3413298.750	96.5388
4		13.270	3427.636	93866.289	2.6548
5		13.270	-4648.466	6300.705	0.1782
Total			263227.728	3535675.983	100.0000