


Sample Description:
Summary:

Product Name: RP10422 (YE-12)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A 70% B 30%
20.0min 100% 0%
25.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10422
Lot No.: PJ00410505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		15.975	618.888	4468.497	0.1914
2		16.232	4199.936	37164.945	1.5915
3		16.847	4247.136	38286.898	1.6395
4		17.123	215286.875	2255291.500	96.5776
Total			224352.835	2335211.840	100.0000