


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient :
 A B
 0.1min 75 25
 25 min 50 50
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10444
 Lot No.: PL00115506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		8.382	5884.120	59660.496	2.0091
2		9.523	1752.254	19044.467	0.6413
3		9.792	3312.411	33481.359	1.1275
4		10.117	204614.859	2857371.000	96.2221
TOTAL			215563.645	2969557.322	100.0000