


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 75 25
 25 min 50 50
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10447
 Lot No.: PL00215506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		8.175	4483.306	37398.004	0.5155
2		8.403	649.130	4891.897	0.0674
3		8.827	5540.984	67793.859	0.9344
4		9.438	386751.563	6894083.000	95.0207
5		10.125	14683.222	251181.844	3.4620
TOTAL			412108.205	7255348.604	100.0000