


Sample Description:

Product Name : RP10456(PEM-11-NH2-1)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10456
 Lot No.: PJ00925504

	A	B
0.01min	25%	75%
25.0min	38%	62%
25.1min	100%	0%
30.0min	STOP	

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.942	2799.172	34792.164	2.0612
2		11.942	1152.532	5704.447	0.3379
3		12.863	393.435	4804.002	0.2846
4		13.317	379.676	3920.645	0.2323
5		13.613	1473.311	15575.490	0.9227
6		13.613	508.496	3376.799	0.2000
7		14.603	106188.320	1617082.250	95.7992
8		15.543	234.550	2736.198	0.1621
Total			113129.494	1687991.994	100.0000