

Sample Description:


Summary:

Product Name: RP10458 (DM-21-NH2)
 Column : 4.6x250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 19% 81%
 25.0min 34% 66%
 25.1min 100% 0%
 30.0min STOP

Cat. No. : RP10458

Lot No. : PX00106502

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		9.443	12774.106	234311.000	0.6781
2		14.332	7678.169	114040.273	0.3300
3		14.578	29710.207	268848.344	0.7780
4		15.240	627852.563	33938252.000	98.2139
Total:			678015.045	34555451.617	100.0000