


Summary:  
 Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient : A B  
           0.1min 65 35  
           25 min 40 60  
           25.01min 0 100  
           30 min 0 100  
 Flow rate : 1.0ml/min

Cat. No.: RP10466  
 Lot No.: PL00415506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		10.492	878.539	18467.596	0.2614
2		11.765	14130.719	239475.375	3.3900
3		12.617	427518.000	6806327.500	96.3486
TOTAL			442527.258	7064270.471	100.0000