


Sample Description:

Product Name : RP10475(RT-5)  
 Column : 4.6\*250mm, Sinochrom ODS-BP 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient  
     A                    B  
 0.01min    10%            90%  
 25.0min    35%            65%  
 25.1min    100%           0%  
 30.0min                    STOP  
 Flow rate : 1.0ml/min  
 Wavelength : 220nm

Cat. No. : RP10475  
 Lot No. : PL00131501



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		15.117	8056.468	80426.898	1.1267
2		16.380	517711.875	6787663.500	95.0845
3		20.117	12680.455	88836.102	1.2445
4		20.407	23087.318	181631.813	2.5444
Total			561536.116	7138558.313	100.0000