


Sample Description:
Summary:

Product Name : RP10488 (PY-24-NH2)
 Column : 4.6*250mm, Hypersil-ODS2 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 25% 75%
 25.0min 50% 50%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP10488
 Lot No.: PX00107503



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		15.105	6200.145	66726.797	0.3554
2		16.172	18924.104	273393.781	1.4563
3		16.388	11831.033	114074.656	0.6077
4		16.852	872160.438	18095022.000	96.3888
5		18.795	7026.449	107267.195	0.5714
6		19.643	8256.532	97963.406	0.5218
7		20.090	1220.496	18500.205	0.0985
Total			925619.196	18772948.041	100.0000