


Sample Description:  
Summary:

Product Name: RP10517(DM-10-NH2)  
Column : 4.6×250mm, Sinochrom ODS-BP 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
0.01min A 30% B 70%  
20.0min A 35% B 65%  
25.1min A 100% B 0%  
30.0min STOP  
Flow rate : 1.0ml/min

Cat. No.: RP10517  
Lot No.: PJ00219505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.358	695.416	9220.398	0.3605
2		14.872	532.618	6940.701	0.2714
3		15.322	482.215	6734.303	0.2633
4		16.153	479.978	13954.794	0.5456
5		16.958	101126.617	2499105.750	97.7085
6		19.245	327.345	4167.502	0.1629
7		19.768	949.973	17592.742	0.6878
Total			104594.162	2557716.191	100.0000