

Sample Description:


Product Name: RP10555[SP-4(1)]
 Column : 4.6×250mm,Sinohrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP10555
 Lot No.: PJ00606506

	A	B
0.01min	3%	97%
25.0min	15%	85%
25.1min	100%	0%
30.0min		

STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.515	1528.019	14617.503	0.4588
2		16.035	214868.141	3029610.000	95.0853
3		17.105	13213.422	141974.156	4.4559
Total			229609.582	3186201.659	100.0000