

Sample Description:

Product Name : RP10566 (ER-8)  
 Column : 4.6\*250mm, Sinochrom ODS-BP 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient


Cat. No.: RP10566

Lot No.: PX00206502

	A	B
0.0 min	5%	95%
25.0 min	15%	85%
25. min	100%	0%
30.0 min	STOP	

Flow rate : 1.0ml/min

Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.107	15998.442	127000.406	0.8553
2		13.397	30846.988	274755.531	1.8505
3		13.647	13683.501	218939.922	1.4745
4		14.168	484139.281	14227269.000	95.8197
Total			544668.213	14847964.859	100.0000