


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 90 10
 25 min 0 100
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10589
 Lot No.: PZ00131501



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		7.183	766730.500	7737427.000	95.2105
2		7.910	5987.729	63448.293	0.7807
3		9.228	24263.998	198504.500	2.4426
4		9.467	9870.162	58895.559	0.7247
5		10.468	4804.415	34400.184	0.4233
6		12.485	823.668	5556.392	0.0684
7		12.672	2528.398	18503.768	0.2277
8		13.055	1209.050	9914.786	0.1220
TOTAL			816217.920	8126650.481	100.0000