


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 90 10
 25 min 0 100
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10589
 Lot No.: PZ00131501



RESULTS

| Peak No. | Peak ID | Ret Time | Height | Area | Conc |
|----------|---------|----------|------------|-------------|----------|
| 1 | | 7.183 | 766730.500 | 7737427.000 | 95.2105 |
| 2 | | 7.910 | 5987.729 | 63448.293 | 0.7807 |
| 3 | | 9.228 | 24263.998 | 198504.500 | 2.4426 |
| 4 | | 9.467 | 9870.162 | 58895.559 | 0.7247 |
| 5 | | 10.468 | 4804.415 | 34400.184 | 0.4233 |
| 6 | | 12.485 | 823.668 | 5556.392 | 0.0684 |
| 7 | | 12.672 | 2528.398 | 18503.768 | 0.2277 |
| 8 | | 13.055 | 1209.050 | 9914.786 | 0.1220 |
| TOTAL | | | 816217.920 | 8126650.481 | 100.0000 |