


Sample Description:

Product Name : RP10689(PYG10-NH2)
Column : 4.6*250mm, Kromasil C18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A 19% B 81%
25.0min A 30% B 70%
25.1min A 100% B 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10689
Lot No.: PJ00529504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.663	8470.688	94766.047	0.6137
2		14.407	16173.586	162909.203	1.0551
3		15.018	981105.875	14934282.000	96.7197
4		15.478	16061.046	174018.203	1.1270
5		16.980	4924.229	74812.195	0.4845
Total			1026735.423	15440787.648	100.0000