


Sample Description:

Structure : RP10692  
 Column : 4.6×250mm, Kromasil C18-5  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient :  
           A          B  
 0.01min 20%      80%  
 25min   30%      70%  
 25.1min 100%    0%  
 30min          STOP

Cat. No.: RP10692

Lot No.: PY00117503

Flow rate : 1.0 mL/min  
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.988	917789.438	22435856.000	95.3102
2		14.363	66286.648	1103980.250	4.6898
Total			984076.086	23539836.250	100.0000