


Sample Description:  
Summary:

Product Name: RP10701(SS-15)  
Column : 4.6×250mm, Kromasil C18 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
                  A                  B  
0.01min          10%              90%  
25.0min          30%              70%  
25.1min          100%            0%  
30.0min                          STOP  
Flow rate : 1.0ml/min  
Wavelength : 220nm

Cat. No.: RP10701  
Lot No.: PL00116503



Results

| Peak No. | Peak ID | Ret Time | Height     | Area         | Conc.    |
|----------|---------|----------|------------|--------------|----------|
| 1        |         | 15.775   | 6338.446   | 90088.648    | 0.6455   |
| 2        |         | 17.490   | 2667.435   | 26869.402    | 0.1925   |
| 3        |         | 17.758   | 23089.346  | 233902.844   | 1.6759   |
| 4        |         | 18.322   | 771957.625 | 13527486.000 | 96.9234  |
| 5        |         | 19.305   | 4193.886   | 66798.102    | 0.4786   |
| 6        |         | 20.800   | 1344.725   | 11733.497    | 0.0841   |
| Total    |         |          | 809591.462 | 13956878.493 | 100.0000 |