


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 65 35
 25 min 40 60
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10702
 Lot No.: PL00815506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		9.410	34.302	10202.634	0.1918
2		9.648	5221.303	61015.465	1.1472
3		10.517	222438.250	5059025.000	95.1208
4		10.517	10053.996	177098.813	3.3298
5		13.420	775.042	11184.496	0.2103
TOTAL			238522.892	5318526.407	100.0000