


Sample Description:

Summary:

Product Name: RP10770(AT-11)
 Column : 4.6×250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 0.01min A B
 25.0min 32% 68%
 25.1min 37% 63%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP10770

Lot No.: PJ00421504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		17.173	2197.771	36778.148	0.6567
2		18.200	193892.031	5380090.500	96.0672
3		18.200	9147.712	183474.781	3.2761
Total			205237.514	5600343.430	100.0000