


## Summary:

Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient :           A           B  
           0.1min 70           30  
           25 min 45           55  
           25.01min 0           100  
           30 min 0           100

Flow rate : 1.0ml/min

Cat. No.: RP10925

Lot No.: PL00408503



## RESULTS

| Peak No. | Peak ID | Ret Time | Height     | Area        | Conc     |
|----------|---------|----------|------------|-------------|----------|
| 1        |         | 4.248    | 7428.595   | 66537.039   | 1.7843   |
| 2        |         | 7.490    | 2998.060   | 34473.176   | 0.9245   |
| 3        |         | 8.752    | 308691.656 | 3585796.250 | 96.1599  |
| 4        |         | 12.103   | 375.676    | 3645.846    | 0.0978   |
| 5        |         | 14.242   | 1368.081   | 38542.125   | 1.0336   |
| TOTAL    |         |          | 320862.068 | 3728994.436 | 100.0000 |