


Summary:

Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 75 25
 25 min 50 50
 25.01min 0 100
 30 min 0 100

Flow rate : 1.0ml/min

Cat. No.: RP10926
 Lot No.: PL00101503



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		11.260	817223.875	16900062.000	99.8057
2		13.952	1262.683	23594.705	0.1393
3		17.370	461.502	9301.512	0.0549
TOTAL			818948.060	16932958.217	100.0000