


Sample Description:

Product Name : RP11099 (CP-32-2(0))
 Column : 4.6*250mm, Hypersil-ODS2 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP11099
 Lot No.: PJ00225505

	A	B
0.01min	30%	70%
25.0min	50%	50%
25.1min	100%	0%
30.0min	STOP	

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.853	5727.867	271706.969	1.8236
2		16.425	13450.106	213614.969	1.4337
3		16.815	500402.469	14414115.000	96.7427
Total			519580.442	14899436.938	100.0000