


Summary:

Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient : A B  
           0.1min 75 25  
           25 min 50 50  
           25.01min 0 100  
           30 min 0 100  
 Flow rate : 1.0ml/min

Cat. No.: RP11116  
 Lot No.: PX00315506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		10.655	693.858	6039.710	0.2815
2		10.905	2723.763	24851.127	1.1584
3		11.548	190844.688	2038760.125	95.0363
4		11.548	10612.141	75592.813	3.5237
TOTAL			204874.449	2145243.775	100.0000