

Sample Description:

Summary:

Product Name: AF-30(O)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient


	A	B
0.01min	21%	79%
25.0min	27%	73%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP11191

Lot No.: PJ00201506

Chromatogram (AF-30-2-3-F-050512.mdy)



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.167	1256.500	16189.900	0.0702
2		15.097	1152.850	12384.999	0.0537
3		15.833	3133.455	52231.449	0.2263
4		17.105	10986.392	161558.578	0.7001
5		18.002	248873.703	22586124.000	97.8691
6		18.002	4349.161	249410.406	1.0807
Total			269752.060	23077899.333	100.0000