


Sample Description:

Product Name : RP11221(QF-13)  
 Column : 4.6\*250mm, Sinochrom ODS-BP 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient

Cat. No.: RP11221  
 Lot No.: PJ01129504

	A	B
0.01min	19%	89%
25.0min	32%	68%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min  
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		18.223	15896.290	367842.000	1.2575
2		18.568	661409.000	28349482.000	96.9184
3		20.143	28891.066	434117.875	1.4841
4		21.323	5458.080	99446.797	0.3400
<b>Total</b>			711654.437	29250888.672	100.0000