


Summary:  
 Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient :           A           B  
           0.1min 85           15  
           25 min 60           40  
           25.01min 0           100  
           30 min 0           100  
 Flow rate : 1.0ml/min

Cat. No.: RP11252  
 Lot No.: PL00312507



### RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		4.945	2016.194	15866.002	0.3064
2		5.523	8587.260	69976.648	1.3515
3		5.802	3634.167	29323.996	0.5663
4		6.147	3614.255	30984.279	0.5984
5		6.497	343918.250	5031653.000	97.1774
TOTAL			361770.126	5177803.926	100.0000