


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 90 10
 25 min 0 100
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP11267
 Lot No.: PX00128502



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		9.605	10278.637	114784.023	1.1851
2		10.343	2043.734	17716.988	0.1829
3		10.565	4311.467	43962.754	0.4539
4		10.870	501763.906	9508920.000	98.1780
TOTAL			518397.745	9685383.766	100.0000