


Sample Description:
Summary:

Product Name : RP11291(SY-40-NH2(0))
 Column : 4.6*250mm, Hypersil-ODS2 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient :
 A B
 0.01min 25% 75%
 25.0min 48% 52%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP11291
 Lot No.: PJ00209510



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		15.902	14031.104	220775.625	2.5008
2		16.523	356265.219	8402682.000	95.1818
3		16.523	19243.592	204579.672	2.3174
Total			389539.915	8828037.297	100.0000