


Summary:

Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 80 20
 25 min 55 45
 25.01min 0 100
 30 min 0 100

Flow rate : 1.0ml/min

Cat. No.: RP11301

Lot No.: PL00204503



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		10.423	12567.307	179523.359	2.6238
2		10.803	477726.688	6613508.000	96.6582
3		11.717	1542.045	15744.385	0.2301
4		12.783	1298.496	33385.898	0.4879
TOTAL			493134.535	6842161.643	100.0000