


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 80 20
 25 min 55 45
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP11302
 Lot No.: PL00615506



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		10.192	13989.302	141727.281	3.4656
2		10.640	316729.250	3939828.500	96.3383
3		12.697	467.710	8021.144	0.1961
TOTAL			331186.262	4089576.925	100.0000